

- For preventing of aerosol contamination of pipette and sample
- For protection of inner parts of the pipette
- Optimized for use with Finnpipettes
- Finntip Filter Micro up to Finntip Filter 1000 µl are guaranteed free of human DNA, DNase and RNase
- Ideal for PCR work
- Irradiated tip racks to guarantee sterility

Finntip Filter




Tip material: virgin polypropylene

- Excellent chemical resistance
- The straight even sides and hydrophobic surface allow complete dispensing of the sample

Filter material: virgin polyethylene

- Non self-sealing filter with hydrophobic nature
- Average pore size 18-40 µm
- Prevents aerosol passage and thus contamination of the pipette and the sample
- Does not contain any sealing additives



Finntip Filter 1000 µl sterile

| Cat.no. | Description | Volume | Quantity |
|-----------|------------------------------------|-------------|-------------|
| 94052 000 | Finntip Filter Micro 10 µl sterile | 0.2-10 µl | 10x96/rack |
| 94052 020 | Finntip Filter Micro 20 µl sterile | 0.5-20 µl | 10x384/rack |
| 94052 060 | Finntip Filter Micro 50 µl sterile | 0.2-50 µl | 10x384/rack |
| 94052 100 | Finntip Filter 10 µl sterile | 0.5-10 µl | 10x96/rack |
| 94052 150 | Finntip Filter 20 µl sterile | 0.5-20 µl | 10x96/rack |
| 94052 160 | Finntip Filter 30 µl sterile | 0.3-30 µl | 10x96/rack |
| 94052 200 | Finntip Filter 100 µl sterile | 0.5-100 µl | 10x96/rack |
| 94052 310 | Finntip Filter 100 µl Ext sterile | 5-100 µl | 10x96/rack |
| 94052 300 | Finntip Filter 200 µl sterile | 0.5-200 µl | 10x96/rack |
| 94052 320 | Finntip Filter 200 µl Ext sterile | 5-200 µl | 10x96/rack |
| 94052 350 | Finntip Filter 300 µl sterile | 5-300 µl | 10x96/rack |
| 94052 410 | Finntip Filter 1000 µl sterile | 100-1000 µl | 10x96/rack |
| 94052 430 | Finntip Filter 1000 µl Ext sterile | 100-1000 µl | 5x96/rack |
| 94052 550 | Finntip Filter 5 ml sterile | 1-5 ml | 5x54/rack |
| 94052 600 | Finntip Filter 10 ml sterile | 2-10 ml | 5x24/rack |

North America: USA/Canada +1 866 984 3766

Europe: Austria +43 1 801 40 0, Belgium +32 2 482 30 30, Finland +358 9 329 100, France +33 2 2803 2000, Germany national toll free 08001-536 376, Germany international +49 6184 90 6940, Italy +39 02 95059 1, Netherlands +31 76 571 4440, Russia/CIS +7 095 225 11 15, Spain/Portugal +34 93 223 3154, Switzerland +41 44 454 12 12, UK/Ireland +44 870 609 9203

Asia: China +86 21 6865 4588 or +86 10 5850 3588, India +91 22 5542 9494, Japan +81 45 453 9220, Other Asian countries +852 2885 4613

Countries not listed: +49 6184 90 6940 or +33 2 2803 2000

www.thermofisher.com